

Introducing TRIOS MOVE Pro



3shape 



Engagement

Ergonomics

Mobility

Meet the new TRIOS MOVE Pro | Designed for Mobility. Powered for Performance.



Always on – built-in battery provides up to 30 minutes of runtime when unplugged*



Improved PC delivers 1.6x more CPU power and 1.5x better GPU performance than MOVE+, ensuring seamless support for TRIOS 6 and its advanced new features



Keep everything you need for a scan in one place – the storage tray holds scanner tips, protection tip, body sleeves, and batteries for convenient access and organized workflow.



Cable hook for easier moveability, keeping cables neatly secured and tangle-free

Meet the new TRIOS MOVE Pro | Engineered for comfort. Designed to enhance engagement.



See it all, upgraded 18.5" touchscreen provides a 40% larger viewing area compared to MOVE+



An effective tool for engaging patients, advancing case acceptance



Adjustable monitor arm to help ensure an optimal, comfortable working position



Sleek design that elevates the look and feel of the clinic

MOVE Pro | Featured Highlights

IMPROVED

18.5" Full HD Touch Display

Experience easier navigation and user-friendly interactions with a large 18.5" FHD touch display.

NEW


Cable Hook

Allowing for easy cable management and a better visual experience when moving the device.

IMPROVED

PC hardware, designed for scanning

Specifically designed for use with our TRIOS 5 and TRIOS 6 scanners to ensure a fast and smooth scan experience. Two USB-C ports.



SAME ERGONOMIC DESIGN

Adjustable Arm

Position the screen directly in front of the patient for greater involvement in the scanning process and treatment planning. The fully adjustable arm enhances ergonomics by allowing the monitor to be placed in any position, enabling the dentist to comfortably sit, stand, or adjust as needed for optimal workflow.

NEW

Storage Tray

Keep everything you need in one place with convenient detachable storage for autoclavable tips or Ready Tips, body sleeves, and batteries. This ensures an efficient and uninterrupted workflow, with all scanning essentials easily accessible.

NEW

Integrated Rechargeable Battery*

Enjoy up to 30 minutes of full mobility between treatment rooms with a built-in battery buffer, ensuring seamless movement between operatories and enhancing workflow flexibility (1).

1) Not during scanning. Scanning will reduce battery runtime.

Portfolio I Technical Overview



MOVE+



MOVE Pro

Features			
Cart Features	Dimensions, H x W x L (cm)	Standing position: 150cm x 50cm x 60cm	Standing position: 150cm x 50cm x 60cm
	Total weight (kg)	13	19
	Display (size)	15.6"	18.5"
	Monitor resolution	Full HD Touch Display	Full HD Touch Display
	Battery	-	✓ 30 minutes
	Speaker	✓	✓ Improved
	Adjustable Monitor Arm	✓	✓
	Scanner Holder	✓	✓
	Storage tray for scanner accessories	-	✓
	Cable Hook	-	✓
PC Specs	Stylus Pen	✓	-
	Processor	Intel® Core™ i5-12600T	Intel i7-14701E
	Graphics	Nvidia T1000 8GB GDDR6	AMD Radeon Pro W7500
	RAM	32GB	32GB
	Storage	1TB SSD	1TB SSD
	Operating System	Windows® 11	Windows® 11
	Wi-Fi	✓	✓ WiFi 5, WPA3
Compatibility	I/O ports	1x USB-A, 1x LEMO	2x USB-C
	TRIOS 6	-	✓
	TRIOS 5	✓	✓
	TRIOS 4	✓	-
	TRIOS 3	✓	-
	TRIOS Core	-	-

MOVE Pro | What's in the box?



Hex key



Hex screw



Hex screwdriver



Scanner holder, incl. two screws



Cloth bag for tools



Storage Tray












Power cable




S&S Guide + additional instructions





2025 Portfolio Overview | MOVE Pro for wireless only

System compatibility		TRIOS 6	TRIOS 5	TRIOS 3 (Wired)	TRIOS Core
					
	MOVE Pro*	✓	✓	✗	✗
	MOVE+	✗	✓	✓	✗
	High-end laptop	✓	✓	✓	✓
	Low-end laptop	✗	✗	✓	✓
	TRIOS Desktop PC	✓	✓	✓	✓

*Note: TRIOS 3 WL and TRIOS 4 are not compatible with TRIOS MOVE Pro

Market availability and expected registration dates



		Country/market	Expected registration
Sept launch		EU	22-08-2025
		Puerto Rico	04-09-2025
		USA	04-09-2025
		Switzerland	30-09-2025
		Great Britain	10-09-2025
		Canada	22-09-2025
Oct launch		Australia	End of October 2025
		Israel	29-10-2025
		Qatar	29-10-2025
		South Korea	End December/beginning Jan
		New Zealand	15-01-2026
		United Arab Emirates	04-02-2026
		Turkey	16-04-2026
		Indonesia	03-08-2026
		Thailand	08-10-2026
		Singapore	30-10-2026
		Kuwait	14-09-2027
		Saudi Arabia	24-04-2026
		Ukraine	10-08-2026

Pricing & TRIOS Care | TRIOS MOVE Pro

TRIOS MOVE Pro

Reseller	5.900 EUR	6.800 USD
End-User	7.000 EUR	8.000 USD

Expected shipment: September 15th 2025



TRIOS MOVE Pro is covered by TRIOS Care extended warranty.



The holder is covered by TRIOS Care extended warranty.



The storage tray is covered by a standard 1-year warranty.



Let's change dentistry together